

---

## PEG Product Information

### PEG-di-Fluorescein

---

Product:	<b><math>\alpha,\omega</math>-Di-fluorescein poly-ethylene glycol</b>
Article Number:	PDF20-1 (100 mg); PDF20-5 (500 mg)
Description:	Aqueous solution of $\alpha,\omega$ -Di-fluorescein poly-ethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG <sub>MW</sub> :	20000 Da
Functional Group:	Fluorescein isothiocyanate conjugate
Excitation:	<b>502 nm</b>
Emission:	<b>525 nm</b>
Storage Buffer:	ddH <sub>2</sub> O
Autoclaved:	Yes
Storage:	4 – 8 °C. <b>PROTECTED FROM LIGHT!</b>
Expiry date:	Two years after production date

---

