
PEG Product Information

PEG-di-Fluorescein

Product:	α,ω-Di-fluorescein poly-ethylene glycol
Article Number:	PDF8-1 (100 mg); PDF8-5 (500 mg)
Description:	Aqueous solution of α,ω -Di-fluorescein poly-ethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG _{MW} :	8000 Da
Functional Group:	Fluorescein isothiocyanate conjugate
Excitation:	502 nm
Emission:	525 nm
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C. PROTECTED FROM LIGHT!
Expiry date:	Two years after production date

