
PEG Product Information

PEG-mono-Phosphate

Product:	α-Methoxy-ω-phosphate polyethylene glycol
Article Number:	PMP5-1 (1 g); PMP5-5 (5 g)
Description:	Aqueous solution of α -Methoxy- ω -phosphate polyethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG _{MW} :	5000 Da
Functional Group:	Phosphate
Storage Buffer:	ddH ₂ O
Autoclaved:	Yes
Storage:	4 – 8 °C.
Expiry date:	Two years after production date

