

---

## PEG Product Information

### PEG-mono-Rhodamine

---

Product:	<b><math>\alpha</math>-Methoxy-<math>\omega</math>-rhodamine polyethylene glycol</b>
Article Number:	PMR5-1 (100 mg); PMR5-5 (500 mg)
Description:	Aqueous solution of $\alpha$ -Methoxy- $\omega$ -rhodamine polyethylene glycol
Weight of Volume:	mg/ml
Lot:	
Production Date:	
PEG <sub>MW</sub> :	5000 Da
Functional Group:	Rhodamine B isothiocyanate conjugate
Excitation:	<b>540 nm</b>
Emission:	<b>625 nm</b>
Storage Buffer:	ddH <sub>2</sub> O
Autoclaved:	Yes
Storage:	4 – 8 °C. <b>PROTECTED FROM LIGHT!</b>
Expiry date:	Two years after production date

---

