

Product Information – SiCORE-Goat anti-mouse IgG

Product:	SiCORE-Goat anti-mouse IgG (whole molecule)
Article Number:	6510-1 (1 ml); 6510-5 (5 ml)
Description:	Aqueous dispersion of silica particles
Application:	For binding of secondary antibodies
Weight of Volume:	10 mg/ml
Lot:	
Production Date:	
Matrix:	Silica, high-porous
Size (hydrodynamic diameter):	1.0 μm
Number of Particles:	$\sim 1.8 \times 10^{12}/\text{g}$
Surface Area:	$\sim 100 \text{ m}^2/\text{g}$
Density:	$\sim 2.25 \text{ g}/\text{cm}^3$
Functional Group:	Goat anti-mouse IgG (whole molecule)
Autoclaved:	No
Storage Buffer:	PBS, 0.05 % sodium azide
Storage:	At 4 – 8 °C. Do not freeze!
Expiry date:	Six months after production date.
Note:	For complete resuspension vortex thoroughly!



— o —

Goat anti-mouse IgG